

Cidox® Application Guidance

Laboratory

Cidox Range	Laboratory Equipment				Laboratory Surfaces			
	Number of tablets	Doses	Litres of water	ppm	Number of tablets	Doses	Litres of water	ppm
Cidox 100	1	1	1.3	75	1	1	1	100
Cidox 1000	1	1	13.3	75	1	1	10	100
Cidox 2000	1	1	27	75	1	1	20	100



Use Cidox according to your organisation's policies and procedures. Store in dry and cool place (25°C or below).

***Guidance only:** Individual site circumstances may vary according to levels of background contamination, temperature, contact time etc.
Use biocides safely. Always read the label and product information before use.